

RESEARCH ONE-PAGER

SS-31

Elamipretide · mitochondria-targeted research

10 mg and 50 mg lyophilized vials

OVERVIEW

SS-31 (elamipretide) is built to sit in the mitochondrial membrane. It is studied for protecting the cell's engines under oxidative stress — energy, fatigue, and tissue that runs hot.

IN PLAIN LANGUAGE

Think of SS-31 as a seatbelt for mitochondria. People look at it when the research question is tired, stressed engines: heart, muscle, nerve, and aging-biology models.

WHAT RESEARCHERS LOOK AT

- Mitochondrial-membrane stability under stress
- Oxidative-stress and energy-failure models
- Fatigue and high-demand tissue research
- A targeted alternative to simply adding more NAD+
- Heart, muscle, and nerve mitochondrial literature

LABORATORY RECONSTITUTION

Reconstitute aseptically with bacteriostatic water. Invert gently until dissolved; do not shake vigorously. A common laboratory starting point is 2 mL bacteriostatic water per 10 mg of peptide (5 mg/mL). Match the volume to the vial on the bench and confirm against the laboratory SOP. Aliquot if the design requires repeated sampling. Offered as 10 mg and 50 mg lyophilized vials.

STORAGE

Lyophilized: refrigerate, protected from light and moisture. After reconstitution: refrigerate and follow laboratory SOP (commonly 14–30 days with bacteriostatic water). Avoid repeated freeze–thaw cycles.

RESEARCH USE ONLY

The language above is a plain-English reading of published research interests. It is not medical advice and not a claim that this product treats, prevents, or diagnoses any condition. This material is not a drug or cosmetic, and it is not intended for human or veterinary use. HarmonyX supplies it to qualified research purchasers for laboratory work. Certificates of analysis are available on request.

Questions on this vial? Paul and Thabby review every request.

research@harmony-x.com · 360-910-6274